

See discussions, stats, and author profiles for this publication at: <https://www.researchgate.net/publication/322580148>

The in vitro release of cytokines and growth factors from fibrin membranes produced through horizontal centrifugation: Characterization of L-PRF-like fibrin membranes

Article in *Journal of Biomedical Materials Research Part A* · January 2018

DOI: 10.1002/jbm.a.36346

CITATIONS

5

READS

1,070

6 authors, including:



Emanuelle Stellet Lourenço

Universidade Federal Fluminense

7 PUBLICATIONS 29 CITATIONS

SEE PROFILE



Carlos Fernando de Almeida Barros Mourão

55 PUBLICATIONS 182 CITATIONS

SEE PROFILE



Paulo Emilio Correa Leite

National Institute of Metrology, Quality and Technology (Inmetro)

33 PUBLICATIONS 532 CITATIONS

SEE PROFILE



Jose Granjeiro

National Institute of Metrology, Quality and Technology (Inmetro)

387 PUBLICATIONS 6,854 CITATIONS

SEE PROFILE

Some of the authors of this publication are also working on these related projects:



Streptococcus pneumoniae [View project](#)



Nanocrystalline Hydroxyapatite Coatings [View project](#)

The *in vitro* release of cytokines and growth factors from fibrin membranes produced through horizontal centrifugation

Emanuelle Stellet Lourenço,¹ Carlos Fernando de Almeida Barros Mourão,¹ Paulo Emílio Corrêa Leite,² José Mauro Granjeiro,^{3,4} Mônica Diuana Calasans-Maia,⁴ Gutemberg Gomes Alves^{5,6}

¹Post-Graduation Program in Dentistry, Fluminense Federal University, Niterói, Rio de Janeiro, Brazil

²Laboratory of Ultrastructure and Cellular Biology Hertha Meyer, Institute of Biophysics Carlos Chagas Filho, Federal University of Rio de Janeiro, Rio de Janeiro, Brazil

³Nacional Institute of Metrology, Quality and Technology, Duque de Caxias, Rio de Janeiro, Brazil

⁴Dentistry School, Fluminense Federal University, Niterói, Rio de Janeiro, Brazil

⁵Department of Molecular and Cell Biology, Institute of Biology, Fluminense Federal University, Niterói, Rio de Janeiro, Brazil

⁶Clinical Research Unit, Antonio Pedro Hospital, Fluminense Federal University, Niterói, Brazil

Received 27 September 2017; revised 17 December 2017; accepted 16 January 2018

Published online 00 Month 2018 in Wiley Online Library (wileyonlinelibrary.com). DOI: 10.1002/jbm.a.36346

Abstract: Platelet-rich fibrin membranes are biomaterials widely used for therapeutic purposes, and canonically produced through the processing of peripheral blood with fixed-angle rotor centrifuges. In this work, we evaluate the *in vitro* stability and release of cytokines and growth factors when these biomaterials are produced with a horizontal swing-out clinical centrifuge. Membranes produced from the blood of 14 donors were morphologically evaluated by scanning electron microscopy and fluorescence microscopy, and their stability was assessed by photographic recording after incubation in culture medium for up to 28 days. The release of 27 cytokines and growth factors was monitored for three weeks through a multiparametric immunoassay. The fibrin membranes presented complex three-dimensional structure with a high density of nucleated cells. A large release of growth factors [platelet derived growth factor, fibroblastic growth factor

(bFGF), and vascular endothelial growth factor] was detected in the first 24 h, followed by time-dependent decay, maintaining significant concentrations after three weeks. Both anti-inflammatory and pro-inflammatory cytokines presented different release peaks, maintaining high rates of elution for up to 21 days. Chemokines of relevance in tissue repair [RANTES, granulocyte colony-stimulating factor (G-CSF)] were also produced in large quantities throughout the experimental period. The present results demonstrate that blood-derived fibrin membranes with high structural stability and cell content can be generated by horizontal centrifugation, being able of a prolonged production/release of growth factors and pro- and anti-inflammatory cytokines. © 2018 Wiley Periodicals, Inc. J Biomed Mater Res Part A: 00B:000–000, 2018.

Key Words: wound healing, fibrin, growth factors, cytokines

How to cite this article: Lourenço ES, Mourão CFAB, Leite PEC, Granjeiro JM, Calasans-Maia MD, Alves GG. 2018. The *in vitro* release of cytokines and growth factors from fibrin membranes produced through horizontal centrifugation. J Biomed Mater Res Part A 2018:00A:000–000.

INTRODUCTION

Due to their high-regenerative potential, platelet concentrates have been widely used in clinical and surgical procedures that require tissue repair,^{1,2} including treatments for chronic wounds, osteonecrosis, infra-osseal defects, and orthopedic procedures.^{3–6} Among those materials, leukocyte- and platelet-rich fibrin (L-PRF) membranes belong to the second generation of platelet aggregates.³ Favored by a rather simple protocol, its production is characterized by a single centrifugation of peripheral blood samples, circumventing the addition of anticoagulants of blood activators to trigger fibrin polymerization.⁷ The confinement of platelets, leukocytes and growth factors in a high-density fibrin network allows it to be handled like a real solid material and

contributes to structural stability for extended periods of time.⁸ This feature turns PRF into a possible scaffold for tissue engineering, with high clinical applicability, simplified protocol, and low cost.⁹

The interesting clinical results of platelet-rich fibrin are directly related to the controlled release of platelet-derived growth factors¹⁰ such as vascular endothelial growth factor (VEGF), platelet derived growth factor (PDGF), fibroblast growth factor (FGF), and transforming growth factor. These factors promote regeneration through processes such as angiogenesis and stimulate the migration and proliferation of osteoblasts, fibroblasts and mesenchymal cells.¹¹ In addition to growth factors, several pro- and anti-inflammatory cytokines can also be produced by leukocytes present in the

Correspondence to: G. G. Alves; e-mail: gutopepe@yahoo.com.br

Contract grant sponsor: Coordenação de Aperfeiçoamento de Pessoal de Nível Superior (CAPES), Fundação de Amparo à Pesquisa do Estado do Rio de Janeiro (FAPERJ), and Conselho Nacional de Desenvolvimento Científico e Tecnológico (CNPq)




FIGURE 1. Fibrin membrane after centrifugation process. (A) Content in three layers: hematocrit, fibrin clot and acellular plasma—from the lower end to the top of the tube. (B) Fibrin membrane rich in platelets and leukocytes after being withdrawn from the tube. (C) Dashed lines indicate the points of sectioning for the cellular evaluation of the membranes: #1 (lower fragment, adjacent to the red blood cell fraction), #2 (medial fragment or main middle portion) and #3 (upper fragment, or the distal portion from the red thrombus).

L-PRF membranes. As these mediators can impact both positively and negatively on the clinical outcome of its use in several treatments, they remain a controversial subject in the scientific literature.^{12,13}

While presenting a rather simple production protocol, there is some controversy in the standardization of the manufacture of leukocyte and platelet-rich fibrin products. Initially developed as an open access technique, nowadays the production of L-PRF is strongly related to the only Food and Drug Administration (FDA)-approved Conformité Européenne (CE)-marked system with certified materials, marketed under the name Intra-Spin L-PRF (Intra-Lock, Boca Raton, FL).^{3,13} This system uses a fixed-angle rotor centrifuge to process a 9 mL aliquot of peripheral blood for 10 min. A recent study investigated the vibration interference from different centrifugation patterns in the quality of PRF membranes obtained with three different models of commercially available table-top centrifuges when compared to the Intra-Spin system. The authors described a deterioration of cellular components and the abnormal formation of the fibrin network with the non-canonical fixed angle centrifuges.¹⁴ However, the literature still lacks on reports of the impact of the production of leukocyte- and platelet- rich fibrin derivatives by the use of horizontal centrifugation with swing-out rotor centrifuges, which are very accessible and of extensive use in clinics and hospitals.

In this context, this study aimed to evaluate the morphological characteristics and *in vitro* structural stability of leukocyte- and platelet- rich fibrin membranes produced with horizontal rotor centrifuges, as well as their capacity of releasing growth factors, as well as pro- and anti-inflammatory cytokines, during the first weeks after production.

MATERIALS AND METHODS

The present work was developed according to the principles recommended for experimentation with human beings determined in the Declaration of Helsinki. The research protocol was approved by the Research Ethics Committee HUAP-UFF (CAAE: 22763513.6.0000.5243). The participants were informed about the procedures and objectives of the study and signed a consent form.

Preparation of platelet rich fibrin membranes

Peripheral blood was collected from 14 healthy donors, aged between 20 and 56 years, with no history of anticoagulant medication use. For the production of each membrane, a 9 mL aliquot of blood, collected in tubes without the addition of any substance (DryVacutube, Biocon, Brazil), was immediately centrifuged at 400 g in a horizontal swing-out rotor centrifuge (B-40, RDE, Brazil) for 10 min. After polymerization of the fibrin matrix, the material was collected, and a slight compression was applied with the aid of a sterile gauze, producing a resistant gelatinous structure, identified as the Platelet Rich Fibrin membrane.

Characterization of the platelet rich fibrin membranes

Scanning electron microscopy. Immediately after production, the membranes were laid over a glass surface, measured with a caliper gauge, and sectioned with a scalpel in three portions of equal length. The parts were identified as: #1 (lower fragment), #2 (medial fragment) and #3 (upper fragment) [Fig. 1(C)], corresponding to the previously described #1, #2 and #3 sections for microscopic evaluation of PRF.¹⁵ The Region 1 was adjacent to the red blood cell fraction (RBC), Region 2 was the main middle portion, and Region 3 was the distal portion from the red thrombus.

Each portion was fixed with Karnovsky's solution and postfixed with 0.2M sodium cacodylate solution and 1% osmium tetroxide, and finally dehydrated in alcohol solutions (ranging from 15 to 100%) and hexamethyldisilazane. The materials were metallized with gold and observed at 15 kV with a scanning electron microscope (JEOL JSM-6490 LV, JEOL, Japan).

Fluorescence microscopy. The distribution of nucleated cells along the membranes was evaluated by fluorescence microscopy. Three samples divided into the previously described #1, #2, and #3 sections were fixed by treatment with 4% paraformaldehyde solution for 15 min. Cell nuclei were evidenced after treatment with 1:5000 4',6-diamidino-2-phenylindole (DAPI), diluted in phosphate buffered saline. The samples were observed with a 20× objective in an inverted fluorescence microscope (Axio.Observer A1, Zeiss, Germany). The microscopic quantitative analysis was




FIGURE 2. Micrographies obtained by SEM of the fragments of fibrin membrane rich in platelets and leukocytes. (A) Upper fragment, region #3. (B) Medial fragment, region #2. (C) Lower fragment, or region #1.

the case of the FGFb, VEGF, and PDGF growth factors, the large release found in the first 24 h was followed by a reduction in concentrations in the subsequent weeks, always significantly lower ($p < 0.05$) after 21 days. A similar pattern of release was obtained for inflammatory cytokines such as IL-12, IL-1 β , TNF- α , and the chemokines GM-CSF and IP-10, although with very different rates of increase for each analyte.

The cytokines IL-1RA and IL-4, with potential anti-inflammatory activity, showed release peaks at the end of the first week, followed by a gradual reduction, but

maintaining levels still higher than the first day of elution. A similar process was observed for the proinflammatory cytokines IFN- γ , eotaxin, and IL-6, the latter presenting a tenfold increase in the release at the 21st when compared to the first day of elution. RANTES, MIP-1a, and MCP-1 also had elution peaks within the first week, but with the maintenance of high release rates throughout the experimental time.

DISCUSSION

In the present work, membranes were produced based on the principles diffused by Dohan et al.,¹³ i.e., the activation of the coagulation cascade due to the platelet shock against the walls of glass tubes in the absence of anticoagulants. Other derivate protocols, employing considerable changes in the G-force allowed the development of other platelet aggregate products, such as A-PRF.¹⁵ Furthermore, for Kobayashi et al.,¹⁶ the change in rotor angulation may impact in the platelet deposition during PRF production. According to those authors, as the rotor angle changes to a more horizontal position, the G force will focus on the tube so as to concentrate the platelet layer in a smaller area, without altering the cell distribution in the membrane, that is, maintaining the concentration of platelets mostly in the buffy coat.

In this context, the present study investigated a protocol employing a swing-out horizontal rotor centrifugal device, widely used in hospital settings such as hematology, and usually accessible in clinical and surgical procedures. Applying this methodology, while maintaining strictly the same G-force (400 G) proposed originally for PRF, the membranes produced presented high similarity in their basic constitution when compared to results from PRF studies.¹⁷ The ultrastructure of the fibrin membrane evidenced by SEM is a dense and intricate gauze-like mesh, very characteristic of the PRF elastic and resistant membranes. It also presents similarity to collagen structures formed *in vivo* during the healing process.¹⁸ The absence of alterations in the size of these membranes (Fig. 4) when kept in culture medium for up to 28 days after production provides a preliminary evidence of stability, similar to that of the PRF in an aqueous medium.^{14,19} It is important to notice that further studies are necessary to confirm that this apparent evidence of stability by direct observation reflects the maintenance on the integrity of the fibrin network, employing techniques such as SEM. Such stability is also a desirable feature for scaffolds for tissue engineering,²⁰ supporting cell migration and proliferation for tissue remodeling purposes, as biodegradable barriers for guided tissue regeneration procedures.²¹

It is worth noting that cell entrapment is a central attribute of leukocyte-rich fibrin membranes. The present fluorescence microscopy assay was able to detect a high density of nucleated cells within the structure of the fibrin membrane, quite typical of PRF membranes, with an uneven distribution along its extension, and a higher cell density at the top portion.²² The presence of apparently viable leukocytes, as well as a fibrin network which seems stable for weeks, could contribute to producing a reservoir for




FIGURE 3. Fluorescence microscopy images with DAPI marker. (A) Lower fragment, or region #1 (Buffy coat). (B) Medial fragment #2. (C) Upper fragment, or region #3. (D) Cell/field count in relation to the membrane segmentation. *Significantly different from the other experimental groups ($p < 0.05$).

controlled release of growth factors and cytokines.^{1,13} In fact, the literature reports the retention and release of these substances from one to four weeks after a surgical intervention with PRF in different tissues.¹

The release of growth factors by cells trapped in the fibrin matrix has been pointed out as a vital element in the regenerative process.^{12,23,24} However, there is controversy regarding the impact of the use of different centrifugation systems in the release of bioactive molecules by membranes produced with otherwise very similar protocols.¹⁴ A widespread approach for the assessment of GF release is the *in vitro* analysis of elution, that is, the change in the concentration of these molecules in a liquid medium. With such methodology, diverse studies have identified the elution of the main growth factors, as well as relevant inflammatory cytokines in culture media. This release may start immediately after membrane production, and going as far as weeks of

incubation of PRF.^{23–26} One study²² was able to identify a relationship between the concentration of growth factors and the region of the membranes, where the highest levels were evidenced in the buffy coat area, similar to the most cellularized fragment of the present work (Figs. 2 and 3), and evidencing the role of leukocytes in the production of growth factors.

In the present work, the *in vitro* elution assay shows that the release of several growth factors and cytokines is maintained in concentrations of biological relevance for weeks after the preparation of blood derived membranes through horizontal centrifugation. The evaluation of a broad range of molecules (27 analytes) was possible by the methodological choice of a multiparametric system based on XMap technology (LuminexCorp). This method allows analyzing multiple analytes in a single sample, increasing the reliability of the results concerning comparisons between




FIGURE 4. Monitoring the structural stability of membranes. (A) Sample photographed after 24 h of incubation. (B) Same sample of figure A, photographed with 28 days of incubation. (C) Evolution of the average length of the samples in millimeters over time. No significant statistical difference was found between experimental times ($p > 0.05$).

TABLE I. Concentrations of Released Mediators on Conditioned Medium After 24 h Elution

Analyte	Concentration (pg/mL) (mean ± SD, n = 8)
FGF-b	19.0 ± 3.1
VEGF	3075.5 ± 1341.7
PDGF-BB	22,0495.6 ± 4687.2
GM-CSF	3.8 ± 1.1
G-CSF	9.8 ± 1.1
IL-1β	77.8 ± 49.2
IL-2	<4.5
IL-6	15.5 ± 21.0
IL-8	12,441.6 ± 8074.9
IL-13	6.5 ± 7.8
IL-15	2.8 ± 0.3
IFN-γ	168.6 ± 15.7
TNFα	57.3 ± 17.5
IL-12p70	65.9 ± 3.4
IL-7	<4.5
IL-17A	<7.1
IL-9	<6.3
IL-5	3.5 ± 0.4
IP-10	4957.5 ± 1656.0
IL-10	9.7 ± 0.4
IL-4	1.2 ± 0
IL-1RA	81.6 ± 35.8
RANTES	1.2 ± 1
MIP-1α	27.8 ± 15.1
MIP-1β	594.5 ± 20.6
MCP-1	131.3 ± 38.8
Eotaxin	84.1 ± 20.4

the relative concentrations of the different analytes, and reduces the chance of conflicting results by experimental variability and manipulation errors.

The role of growth factors released by fibrin membranes on tissue regeneration is well known, either through the promotion of angiogenesis, migration, and proliferation of mesenchymal cells and fibroblasts and osteoblasts.¹⁰ In this study, significant levels of released VEGF, bFGF, and PDGF-BB were identified in the first 24 h, similarly to reports on PRF membranes.^{7,19,23} However, methodological differences such as the volume of medium used for elution, detection technique (ELISA, flow luminometry) as well as the parameters for normalization of results (e.g., membrane mass, initial volume of blood samples) prevent an adequate comparison of the absolute concentrations identified from different reports. Nevertheless, a gradual reduction in the release of growth factors over time is confirmed in all of these studies, and in agreement with the present results. It is worth to note that, after 21 days, the released concentrations of molecules such as PDGF correspond to biologically relevant levels identified in previous studies.²⁷

The interpretation of the impact of released molecules by the fibrin membranes and similar becomes more difficult when considering other molecules not identified as growth factors, such as the chemokine RANTES or CCL-5. A comparative study of different platelet derivatives²³ showed that the release of CCL-5 is much higher in PRP than in PRF in

the first 24 h after its production. In fact, the present result for 24 h indicates a low release of this molecule but followed by a significant increase (more than 10-fold) from the 7th to the 21st day. While the possible effects of an increased release of RANTES in the clinical outcome of PRF remain undescribed, there is evidence in the literature of RANTES acting positively in angiogenesis, as well as potential effects of this molecule on the differentiation of osteoblasts during tissue regeneration.²⁸

In the present study, high concentrations of inflammatory cytokines such as IL-1β, IL-6, IL-5, IL-15, TNF-α, and IFN-γ were detected from the first to the 21st day, indicating extensive activation of immune cells in the material. In fact, the presence of immune system cells within the fibrin membrane, and the consequent release of inflammatory mediators has been a subject of controversy in the literature. Although proinflammatory cytokines such as IL-6 may affect the regeneration process of lesions,¹¹ they also play important immunological functions against possible infections and help to clean the perimeter of the wound. According to Dohan et al.,⁷ the leukocyte-rich PRF membrane may represent an immune nodule capable of stimulating defense mechanisms. The high release of proinflammatory cytokines concomitant with growth factors could counterfeit a possible osteolysis⁷ by the simultaneous induction of proliferation of osteoblasts, fibroblasts, and endothelial cells. Such duality may result in a rapid remodeling without the classic signs of the inflammatory picture (heat, flushing, pain, edema, and loss of function), as observed in clinical trials of PRF.²⁹ This may also be related to the long-term release of anti-inflammatory cytokines such as the antagonist of the IL-1 receptor (IL-1RA), which was consistently released in the present study, suggesting the possibility of induction of




FIGURE 5. Heatmap of the variation of analyte concentrations in culture media from 1 to 21 days of incubation of fibrin membranes. The analytes IL-2, IL-7, IL-9, and IL-17 had undetectable concentrations.

a self-limited inflammatory process. Also, *in vitro* studies points to an active role for some proinflammatory cytokines in tissue repair, including evidence that IL-1 β released by PRF membranes may be associated with the migration of mesenchymal stem cells and human endothelial cells.¹⁹ Furthermore, the high concentrations of G-CSF found in this work may also reinforce the potential for promotion of tissue remodeling, since effects on bone regeneration via revascularization and osteogenesis have been described this molecule.³⁰

As recently stated by Dohan Ehrenfest et al.,¹⁴ changes in the protocol of production of L-PRF, using centrifugation systems other than the CE/FDA approved Intra-Spin L-PRF system often result in the biological signature of the resulting material. In this context, the clinical and experimental results from the original L-PRF shall not be extrapolated to the material produced in the present protocol. Nevertheless, taken together, the results of the present work demonstrate that the structural stability of fibrin membranes generated by horizontal centrifugation allows the entrapment of nucleated cells that might contribute to a consistent release of modulatory molecules for weeks after their production. Therefore, further studies are needed to confirm whether this *in vitro* release is parallel to the physiological environment, as well as its actual impact on the clinical outcome of treatments involving these growth factor concentrates. Nevertheless, the present results demonstrate that blood derived fibrin membranes with high structural stability and cell content can be produced by horizontal centrifugation, which are able of a prolonged production/release of growth factors and pro- and anti-inflammatory cytokines.

REFERENCES

- Borie E, Oliví DG, Orsi IA, Garlet K, Weber B, Beltrán V, Fuentes R. Platelet-rich fibrin application in dentistry: A literature review. *Int J Clin Exp Med* 2015;8:7922–7929.
- Anitua E, Pino A, Orive G. Opening new horizons in regenerative dermatology using platelet-based autologous therapies. *Int J Dermatol* 2017;56:247–251.
- Dohan Ehrenfest DM, Andia I, Zumstein MA, Zhang C-Q, Pinto NR, Bielecki T. Classification of platelet concentrates (platelet-rich plasma-PRP, platelet-rich fibrin-PRF) for topical and infiltrative use in orthopedic and sports medicine: Current consensus, clinical implications and perspectives. *Muscles Ligaments Tendons J* 2014;4:3–9.
- Piccin A, Di Pierro AM, Canzian L, Primerano M, Corvetta D, Negri G, Mazzoleni G, Gastl G, Steurer M, Gentilini I, Eisendle K, Fontanella F. Platelet gel: A new therapeutic tool with great potential. *Blood Transfus* 2016;1–8.
- Joshi Jubert N, Rodríguez L, Reverté-Vinaixa MM, Navarro A. Platelet-rich plasma injections for advanced knee osteoarthritis: A prospective, randomized, double-blinded clinical trial. *Orthop J Sport Med* 2017;5:232596711668938. 2325967116689386.
- Kim J-W, Kim S-J, Kim M-R. Leucocyte-rich and platelet-rich fibrin for the treatment of bisphosphonate-related osteonecrosis of the jaw: A prospective feasibility study. *Br J Oral Maxillofac Surg* 2014;52:854.
- Dohan DM, Choukroun J, Diss A, Dohan SL, Dohan AJ, Mouhyi J, Gogly B. Platelet-rich fibrin (PRF): A second-generation platelet concentrate. Part III: Leucocyte activation: A new feature for platelet concentrates? *Oral Surg Oral Med Oral Pathol Oral Radiol Endod* 2006;101:e51–e55.
- Ahmed TAE, Dare EV, Hincke M. Fibrin: A versatile scaffold for tissue engineering applications. *Tissue Eng. Part B Rev* 2008;14:199–215.
- Gassling V, Douglas T, Warnke PH, Agil Y, Wiltfang J, Becker ST. Platelet-rich fibrin membranes as scaffolds for periosteal tissue engineering. *Clin Oral Implants Res* 2010;21:543–549.
- Dohan Ehrenfest DM, de Peppo GM, Doglioli P, Sammartino G. Slow release of growth factors and thrombospondin-1 in Choukroun's platelet-rich fibrin (PRF): A gold standard to achieve for all surgical platelet concentrates technologies. *Growth Factors* 2009;27:63–69.
- De Pascale MR, Sommese L, Casamassimi A, Napoli C. Platelet derivatives in regenerative medicine: An update. *Transfus Med Rev* 2015;29:52–61.
- Anitua E, Zaldueño M, Troya M, Padilla S, Orive G, da Costa Martins PA. Leukocyte inclusion within a platelet rich plasma-derived fibrin scaffold stimulates a more pro-inflammatory environment and alters fibrin properties. *PLoS One* 2015;10:e0121713.
- Dohan DM, Choukroun J, Diss A, Dohan SL, Dohan AJ, Mouhyi J. Platelet-rich fibrin (PRF): A second-generation platelet concentrate. Part I: Technological concepts and evolution. *Oral Surg Oral Med Oral Pathol Oral Radiol Endod* 2006;101:e37–e44.
- Dohan Ehrenfest DM, Pinto NR, Pereda A, Jiménez P, Corso MD, Kang BS, Nally M, Lanata N, Wang HL, Quirynen M. The impact of the centrifuge characteristics and centrifugation protocols on the cells, growth factors, and fibrin architecture of a leukocyte- and platelet-rich fibrin (L-PRF) clot and membrane. *Platelets* 2017;1–14.
- Kobayashi M, Kawase T, Horimizu M, Okuda K, Wolff LF, Yoshie H. A proposed protocol for the standardized preparation of PRF membranes for clinical use. *Biologicals* 2012;40:323–329.
- Swamy PM. *Laboratory Manual on Biotechnology*, Rastogi Publications; 2008
- Ghanaati S, Booms P, Orlowska A, Kubesch A, Lorenz J, Rutkowski J, Landes C, Sader R, Kirkpatrick C, Choukroun J. Advanced platelet-rich fibrin: A new concept for cell-based tissue engineering by means of inflammatory cells. *J Oral Implantol* 2014;40:679–689.
- Dohan Ehrenfest DM, Del Corso M, Diss A, Mouhyi J, Charrier J-B. Three-dimensional architecture and cell composition of a Choukroun's platelet-rich fibrin clot and membrane. *J Periodontol* 2010;81:546–555.
- Wu C-L, Lee S-S, Tsai C-H, Lu K-H, Zhao J-H, Chang Y-C. Platelet-rich fibrin increases cell attachment, proliferation and collagen-related protein expression of human osteoblasts. *Aust Dent J* 2012;57:207–212.
- Schär MO, Diaz-Romero J, Kohl S, Zumstein MA, Nestic D. Platelet-rich concentrates differentially release growth factors and induce cell migration *in vitro*. *Clin Orthop Relat Res* 2015;473:1635–1643.
- Galav S, Chandrashekar KT, Mishra R, Tripathi V, Agarwal R, Galav A. Comparative evaluation of platelet-rich fibrin and autogenous bone graft for the treatment of infrabony defects in chronic periodontitis: Clinical, radiological, and surgical reentry. *Indian J Dent Res* 2016;27:502–507.
- Bansal M, Kumar A, Puri K, Khatri M, Gupta G, Vij H. Clinical and histologic evaluation of platelet-rich fibrin accelerated epithelization of gingival wound. *J Cutan Aesthet Surg* 2016;9:196–200.
- Nishimoto S, Fujita K, Sotsuka Y, Kinoshita M, Fujiwara T, Kawai K, Kakibuchi M. Growth factor measurement and histological analysis in platelet rich fibrin: A pilot study. *J Maxillofac Oral Surg* 2015;14:907–913.
- Passarelli F, Tia M, D'Esposito V, Pascale MD, Corso MD, Sepulveres R, Liguoro D, Valentino R, Beguinot F, Formisano P, Sammartino G. Growth-promoting action and growth factor release by different platelet derivatives. *Platelets* 2014;25:252–256.
- Kobayashi E, Flückiger L, Fujioka-Kobayashi M, Sawada K, Sculean A, Schaller B, Miron RJ. Comparative release of growth factors from PRP, PRF, and advanced-PRF. *Clin Oral Investig* 2016;20:2353–2360.
- Qiao J, An N, Ouyang X. Quantification of growth factors in different platelet concentrates. *Platelets* 2017;28:774–778.
- DuPont NC, Wang K, Wadhwa PD, Culhane JF, Nelson EL. Validation and comparison of luminex multiplex cytokine analysis kits

- with ELISA: Determinations of a panel of nine cytokines in clinical sample culture supernatants. *J Reprod Immunol* 2005;66:175–191.
28. Colciago A, Celotti F, Casati L, Giancola R, Castano SM, Antonini G, Sacchi MC, Negri-Cesi P. In vitro effects of PDGF isoforms (AA, BB, AB and CC) on migration and proliferation of SaOS-2 osteoblasts and on migration of human osteoblasts. *Int J Biomed Sci* 2009;5:380–389.
 29. Czekanska EM, Ralphs JR, Alini M, Stoddart MJ. Enhancing inflammatory and chemotactic signals to regulate bone regeneration. *Eur Cell Mater* 2014;28:320–334.
 30. Jankovic S, Aleksic Z, Klokkevold P, Lekovic V, Dimitrijevic B, Kenney EB, Camargo P. Use of platelet-rich fibrin membrane following treatment of gingival recession: A randomized clinical trial. *Int J Periodontics Restorative Dent* 2012;32:e41–e50.
 31. Ishida K, Matsumoto T, Sasaki K, Mifune Y, Tei K, Kubo S, Matsushita T, Takayama K, Akisue T, Tabata Y, Kurosaka M, Kuroda R. Bone regeneration properties of granulocyte colony-stimulating factor via neovascularization and osteogenesis. *Tissue Eng Part A* 2010;16:3271–3284.